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Optimum design and efficient utilization of catalysts require a thorough understanding of the surface structure and surface chemistry of the catalytic material. Chemical adsorption (chemisorption) analyses can provide much of the information needed to evaluate catalyst materials in the design and production phases, as well as after a period of use. The chemical adsorption isotherm reveals information about the active surface of a material. In addition, temperature-programmed reaction techniques have emerged as an indispensable companion to chemisorption isotherm analyses in many areas of industry and research.

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
AutoChem II 2920 Chemisorption Analyzer is a fully automated instrument capable of conducting a comprehensive array of highly precise chemical adsorption and temperature-programmed reaction studies. The instrument enables the researcher to obtain valuable information about the physical properties of catalysts, catalyst support, and a variety of other materials as shown above.

実用例 / Application Example

Condition: The sample was pretreated at 200°C for 0.3 h under He (50 mL/min) atmosphere. The reduction was carried out from 80 to 1000°C (10°C/min) under 5 vol% H₂ in Ar with a flow rate of 50 mL/min.

Figure 1 TPR (temperature programmed reduction) spectra of Ni-oxide, Sn-oxide



Figure 2 TCD (Thermal Conductivity Detector) signals of ZSM-5 and ammonia

Pulse Chemisorption Report

Experiment 1: Pt Alumina
Analysis Type: Pulse Chemisorption
Calibration: None
Measured Flow Rate: 49.99 cm³ STP/min
Signal Offset: 0.12525
Signal Inverted: No

Peak Number	Temperature at Maximum (°C)	Quantity Adsorbed (cm ³ STP)	Cumulative Quantity (cm ³ STP)
1	34.6	0.19835	0.19835
2	34.4	0.00591	0.20426
3	34.4	0.00548	0.20974
4	34.5	0.00098	0.21072
5	34.5	0.00000	0.21072
6	34.5	0.00000	0.21072

Pulse Chemisorption Analysis Summary

Element	Percent of Sample Mass (%)	Atomic Weight	Stoichiometry Factor	Atomic Cross-Sectional Area (nm ²)	Density (g/cm ³)
platinum	0.5000	195.080	1.000	0.0800	21.450

Active Loop Volume at 110.8 °C: 0.15223 cm³ STP
Cumulative Quantity: 0.20540 cm³ STP
Metal Dispersion: 35.7551%
Metallic Surface Area: 0.4415 m²/g sample
Metallic Surface Area: 88.3052 m²/g metal
Active Particle Diameter: 3.1675 nm
Cubic Crystallite Size: 2.636 nm

Figure 3 CO pulse chemisorption on 0.5 Pt Almina under various temperature