

電界放出形走査型電子顕微鏡

Field Emission Scanning Electron Microscope (FE-SEM)

基本仕様 / Specifications

日立ハイテック (Hitachi High-Tech) SU6600

- 加速電圧 / Acceleration Voltage : 0.5 ~ 30 kV
- 二次電子分解能 / Secondary electron resolution
 - 高真空モード / High vacuum mode : 1.2 nm (30 kV), 3.0 nm (1 kV)
 - 低真空モード / Low vacuum mode (150 Pa) : 4.5 nm (30 kV)
- 反射電子分解能 / Backscattered electron resolution : 3.0 nm (30 kV, 低真空 (10 Pa))
- エネルギー分散型X線分析装置搭載 / Energy dispersive X-ray spectrometer (EDX) XFlash® 5010, Bruker
- 電子後方散乱回折装置搭載 / Electron backscatter diffractometer (EBSD) C-Nano+, Oxford Instruments



基本原理 / Mechanism





Fig. 1 試料から得られる情報 (Information obtained from specimen)

分解能 (Resolution)	200 nm	0.1 nm	0.5 nm
倍率 (Magnification)	~ ×2000	× 50 ~ 1,500,000	× 10 ~ 1,000,000

Scanning electron microscopy (SEM) enables the precise observation of microscopic surface structures, which is impossible with an optical microscope. Moreover, as it can provide images with deeper focal depths, it enables the observation of three-dimensional images, similar to observing a substance with the naked eye, by enlarging the specimen surface, which has a rough structure.

SEM forms an image by using electrons (secondary electrons in general) that are reflected or generated from the surface of the specimen. Because the intensity of the generated secondary electrons varies depending on the angle of the incident electrons on the specimen surface, subtle variations in the roughness of the surface can be expressed according to the signal intensity.

In addition, elemental information can be obtained by analyzing characteristic X-rays, which are detected using energy-dispersive X-ray spectroscopy (EDX). Furthermore, electron backscatter diffraction (EBSD) can detect the diffraction patterns of electron beams that are specific to the crystal structure and orientation near the sample surface, thereby providing crystallographic information regarding the sample microstructure.

<https://www.jeol.co.jp/science/sem.html>

実用例 / Application Examples




Figure 1. FE-SEM micrograph of the surface of bioactive polymeric bone implant spontaneously covered with flake-like crystallites of hydroxyapatite, which is a main component of bone tissue, in simulated body environment.



Figure 2. FE-SEM micrograph of microcapsules consisted of flake-like crystallites of hydroxyapatite formed in simulated body environment.




Figure 3. FE-SEM and EDX mapping integrated image of cross-section of self-supported thin film consisted of hydroxyapatite (HAp; $\text{Ca}_{10}(\text{PO}_4)_6(\text{OH})_2$) formed in simulated body environment.